




Improving Gradient Accuracy and Precision in HPLC Instrumentation



ENGINEERING YOUR SUCCESS.

Ultra Low Carryover Gradient Valve Miniature Liquid Valve


Reduce instrument size and weight

Reduced Size

An 80% size reduction and 75% decrease in weight vs. competing gradient valves with superior flow capabilities

Increase Design Flexibility

Space occupied by the Ultra Low Carryover Gradient Valve is about 80% less than competing valves



Reduced Dwell Volume

The Ultra Low Carryover Gradient Valve reduced size enables the valve to be placed closer to the injector reducing dwell volume and increasing throughput

Enables Smaller and Lighter Systems

This can reduce footprint and instrument weight

Ultra Low Carryover Gradient Valve Miniature Liquid Valve


Improve fluidic performance & improve change over between solvent channels

The Ultra Low Carryover Gradient Valve reduces internal volume in two ways:


- Reducing the overall size of the valve shortens fluid passageways
- The unique armature design extends into the manifold body further reducing internal volume

Channel volume from seal to convergence point is reduced by 40%


- Improves gradient responsiveness
- This low carryover design reduces crosstalk between channels improving gradient accuracy



Ultra Low Carryover
Gradient 2-Channel



Ultra Low Carryover
Gradient 4-Channel



Competing Valve

Ultra Low Carryover Gradient Valve Miniature Liquid Valve

How the electronics design of the Ultra Low Carryover Gradient Valve improves gradient performance

Patented Solenoid


The patented solenoid design enables extremely low response time of 2 msec or less enabling rapid and precise gradient changes

Accuracy

Superior response time and repeatability results in improved compositional accuracy


Low Response Time

When tested with competing valves, the Ultra Low Carryover Gradient Valve showed a 60% reduction in response time



Response Time and Low Internal Volume

This best-in-class response time in combination with the lower internal volume of the Ultra Low Carryover Gradient Valve, enable the best compositional accuracy and precisions



Ultra Low Carryover Gradient Valve Miniature Liquid Valve

The Ultra Low Carryover Gradient Valve offers unparalleled long life to your system.

Life Rating

With a life rating of 50 million cycles the Ultra Low Carryover Gradient Valve offers five times the life of competing valves.




Replacement Interval

With this remarkable life rating this valve can last the life of the HPLC instrument

Reduce Warranty Cost

This improved reliability will reduce failures within the HPLC system warranty period decreasing warranty costs and downtime



Ultra Low Carryover Gradient Valve Miniature Liquid Valve

The Ultra Low Carryover Gradient Valve provides the best crystallization resistance to your design

Saline Tests


In side by side testing the Ultra Low Carryover Gradient Valve still passed specifications after 30 Million cycles with saline solution*

Outperforms Competition

Competing valves failed before their 10M cycle rating

Reduce Downtime

The ability to pass crystals without impacting valve performance reduces warranty and downtime costs




* Valve was still meeting all specifications when test was ended at 30 million cycles

Ultra Low Carryover Gradient Valve Miniature Liquid Valve

Fluidic System Simplification

The Ultra Low Carryover Gradient Valve can also be used as a solvent selector valve in either 2 or 4 channel configurations.



Three different designs to meet your needs:



**4-Channel
Radial Design**



**4-Channel
Panel Mount**

(Can be installed on a panel separating the liquid connections from the electronics)



2-Channel Design


Ultra Low Carryover Gradient Valve Miniature Liquid Valve

Product Specifications

Physical Properties		Electrical		Performance Characteristics	
Valve Type:		Voltage (VDC ±):	12	24	Leak Rate:
4-Channel & 2-Channel		Power (Watts):	3.0	3.0	0.150 sccm of Air (bubble tight)
Media:		Current (mA):	250	116	Operating Pressure:
Liquid		Resistance (Ohm):	48	207	45 psig (3.1 bar)
Operating Environment/ Media Temperature:		$\Omega \pm 10\% @ 68^\circ\text{F}, 20^\circ\text{C}$ Note: For actuation exceeding 100ms Hit & Hold is required.		Proof Pressure:	
4°C to 50°C (39°F to 122°F)				Response Time:	
Storage Temperature:		Electrical Termination:	Molex Housing: #50-57-9402 Molex Contacts: #16-02-0098		120 psig (8.3 bar)
-20°C to 70°C (-4°F to 158°F)		Wetted Materials*		<2 msec at 2X rated Voltage <10 msec at rated Voltage	
Weight:		Seals:	FFKM or EPDM		
4-Channel Radial Design	2.56oz (72.5g)	Body:	PEEK		Recommended Filtration:
4-Channel Panel Mount	2.45oz (69.5g)	Regulatory:	RoHS directive (2002/95/EC) and REACH EC 1907/2006		16 µm
2-Channel Design	1.19oz (33.7g)	Mounting Options:	1/4-28 Threaded Female Design		Reliability:
Porting:					50 Million Cycles
1/4-28					Flow Rate:
Internal Volume: (seal to common port)					Minimum water flow of 320 mL/min @ 45 psig (3.1 bar)
4-Channel Radial Design	9.4 µL				
4-Channel Panel Mount	12.99 µL				
2-Channel Design	4.05 µL				
Orifice Size:					
.030" (0.76 mm)					

ULC	3	24	FF	3	R	F	-000
Series	Configuration	Voltage	Seal Manifold	Orifice	Mounting	Electrical Connection	Configuration
ULC-	3: 2-Channel	12: 12 VDC 24: 24 VDC	FF: FFKM EP: EPDM	3: 0.030" (0.76mm)	4: 1/4 - 28	F: Latching Connector	-000
ULC-	5: 4-Channel	12: 12 VDC 24: 24 VDC	FF: FFKM EP: EPDM	3: 0.030" (0.76mm)	S: Panel Mount 1/4-28 R: Radial Body 1/4-28	F: Latching Connector	-000
Accessories							
Part Number			Description				
290-006061-005			19.5 in (495.3mm) Wire Extension with Flying Leads				

Parker, your partner in fluidic circuit development



With over 30 years of expertise in integrating fluidic circuits, Parker is in a unique position to assist you with your instrument designs. We are the only company that manufactures liquid valves, liquid pumps, gas valves, and gas pumps. Because we manufacture both pumps and valves, you can rely on our expertise to provide a reliable and cost-effective solution. This expertise helps solve your fluidic needs by providing products or integrating them onto manifolds. If you are looking for a pre-tested solution, the entire subsystem can be delivered as a module.

Clinical Diagnostics and Analytical Instrumentation Expertise

Liquid Valves

Miniaturized valves featuring inert materials for the highest chemical compatibility, long life, low carryover and high pressure.

40
YEARS

Material Science – Sealing Elastomers and Valve Bodies

Elastomers and other wetted path materials developed for improved chemical compatibility, long life and optimized temperature performance.

30
YEARS

Customization for OEM Projects

As your valve and pump engineers, we can optimize valve and pump performance to suit your OEM application.

40
YEARS

Ultra Low Carryover Gradient Valve Miniature Liquid Valve

Serving a broad spectrum of life science, air quality, and process instrumentation OEM fluidic needs



Providing Pressure and Vacuum:
Broad range of diaphragm pumps for Gas and Liquid



Gas Flow Control:
High to Low Flow Proportional Valves



On/Off & Channel Selection Capabilities:
Gas and Liquid Solenoid Valves



High Precision Thermal Flow Control:
Mass Flow Controllers and Meters

In order to provide the best possible solution for your application, please provide the following requirements when contacting Applications Engineering:

- Media, Inlet & Outlet Pressures
- System Supply Voltage
- Ambient Temperature Range
- Minimum Required Flow Rate
- Media and Media Temperature Range
- Carryover requirement and how measured

For more information call +1 603 595 1500 or email ppfinfo@parker.com

Visit www.parker.com/precisionfluidics

Recommendations on application design and material selection are based on available technical data and are offered as suggestions only. Each user should conduct their own tests to determine the suitability for their own use. Parker offers no express or implied warranties concerning the form, fit, or function of a product in any application.