


## Configuration


1" Hosebarb

3/4" Tri-Clamp

1" Tri-Clamp


1 1/2" Tri-Clamp

1" Tri-Clamp (T-Port)

## Dimensions

Cartridge Type	Dimension 'A'	Dimension 'B'
10" (250 mm)	10.30" (262 mm)	13.07" (332 mm)
20" (500 mm)	20.04" (509 mm)	22.79" (579 mm)
30" (750 mm)	29.80" (757 mm)	32.56" (827 mm)

Dimensions shown are typical lengths for 1 1/2" Tri-Clamp.  
Further dimensions available from your local representative.

## English

### Introduction

These guidelines give the correct methods for installing and operating our Capsules. If you have any queries please contact our process filtration specialists who will be pleased to discuss your particular filtration requirements or answer any questions you may have.

### 1. Storage

Store Capsules in a clean and dry environment and avoid placing heavy objects on top of the Capsule box. Capsules should not be stored at temperatures below 5 °C (41 °F) or above 30 °C (86 °F) or in direct sunlight.

#### 1.1 The shelf-life for Capsules is as follows:

Membrane capsules, non-gamma irradiated	3 years
Depth capsules, non-gamma irradiated	5 years
Membrane / depth capsules, gamma irradiated	2 years

### 2. Operating Conditions

We certify that this product complies with the European Council Pressure Equipment Directive (PED) 2014/68/EU Article 3, Paragraph 3 – Sound Engineering Practise (SEP). This product is intended for use with Group 1 & 2 Dangerous and Harmless Liquids and Group 2 Harmless Gases at the following maximum operating conditions :

Liquid:	Maximum 5.5 barg (79.7 psig) @ 25 °C ( 77 °F)
	Maximum 2.8 barg (40.6 psig) @ 60 °C (140 °F)
Gas:	Maximum 4.0 barg (58.0 psig) @ 25 °C ( 77 °F)

In compliance with PED Article 3, Paragraph 3, SEP, this product does not bear the CE mark.

### 3. Installation and Operation

Please ensure that the Capsule you have chosen is compatible with the fluid and system before installation and check that the product and part number etched on the base of the Capsule are as required.

**NB: If the capsule is to be used in a sterile environment please ensure that the capsule is either sterile packed or autoclaved prior to use. Capsules can not be steam sterilized in-situ.**

#### a) Liquid Systems

- Follow the arrows on the Capsule label and ensure that the Capsule is installed in the correct direction and using the appropriate connections.
- Before passing liquid through the capsule ensure that the vent valves are closed and the Capsule is mounted in a secure position. It is important that the Capsule is full of liquid during operation, so ensure that it is mounted so that the air can be vented in a safe manner.
- The Capsule can be vented by turning the vent nut in an anti-clockwise direction. The valve is retained by the nut so do not undo the nut by more than 1 or 2 turns as this may result in wastage of the liquid within the Capsule or injury to the operator.  
**NB: If hazardous fluids are to be passed through the Capsule, the vent valves should be connected to an appropriate drain line.**
- If the Capsule is connected into a pressurized system make sure that transient pressures do not exceed the maximum operating conditions.

#### b) Gas Systems

- Follow the arrows on the capsule label and ensure that the Capsule is installed in the correct direction and using the appropriate connections. For vent applications flow can be in either direction.
- Gas filters are hydrophobic and may not allow the passage of water. It is also possible that condensation can collect in the Capsule and prevent the passage of gas. The Capsule should be mounted so that any water or condensate can be drained from the vent valves. In some cases it may be necessary to remove the Capsule and dry the membrane by passing dry air or nitrogen through it.

### 4. Integrity Testing

Some Capsules can be integrity tested either off-line or in-situ. Please contact your local representative for further information on the most appropriate method and equipment or refer to the appropriate product datasheets.

### 5. Sterilization

Capsules can be autoclaved up to the maximum temperature marked on the product. Please contact your local representative with the details of your sterilization regime so we can advise the maximum number of cycles.

- NB:**
- The vent / drain valves should be left open during the autoclave cycle
  - DO NOT STEAM STERILIZE CAPSULES IN-LINE !

### 6. Disposal

All Capsules should be disposed of in a safe manner and in line with Health and Safety Guidelines.

## Installation and Operating Instructions for Capsules

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